

CTN

RECOVERY TECHNOLOGIES

LedProTM

Multimodal Wellness System

RED/NIR LIGHT EXPOSURE
VIBRATION TRAINING
NOGIER-INSPIRED SOUND FREQUENCIES

USER MANUAL



CTN


26.08.2026

TABLE OF CONTENTS


TECHNICAL SPECIFICATIONS	2
DEVICE SCHEMATICS	3
INTRODUCTION	4
CONTRAINDICATIONS	4
SAFETY PRECAUTIONS	5
UNBOXING INSTRUCTIONS	10
INTRODUCTION TO FUNCTIONS	11
RED LIGHT AND NEAR-INFRARED	11
VIBRATION TRAINING	11
NOGIER-INSPIRED SOUND FREQUENCIES	11
INSTRUCTIONS FOR USE	12
STARTUP PROCEDURE, STARTING AND STOPPING THE SESSION	12
SESSION GUIDELINES	13
SESSION PRESETS	13
MANUAL SESSION SETTINGS	14
MAINTENANCE	15
WARRANTY INFORMATION	15

TECHNICAL SPECIFICATIONS

Product Name:	LedPro™ Vertical
Display:	24" full HD monitor / LCD touchscreen
Power Consumption:	3200W
Panel Power:	1200 W per panel
Power Supply (AC):	230 V / 50 Hz, 16 A fuse required
LED Amount:	840 per panel (420 red ; 420 near-infrared)
Cooling System:	Dual AC radial fans
Main Light Wavelengths:	660 nm (±10) red + 850nm (±10) NIR
Total Irradiance:	160 mW/cm² @ 15-30 cm distance, per panel
Vibration Plate Power:	380 W
Vibration Plate Frequency Range:	1-40 Hz
Vibration Plate Amplitude:	0-5 mm @ 57,5 kg load
Sound System:	JOY Soundmate 3
Max User Weight:	240 kg
Device Weight:	310 kg
Device Dimensions:	860 x 1150 x 2240 mm (L*W*H)
Transport box dimensions:	2377 x 998 x 1426 mm (L*W*H)
Session Modes:	BODY TONING, SKIN TONING, INFLAMMATION REDUCTION, MOBILITY, STRESS RELIEF



DEVICE SCHEMATICS



INTRODUCTION

Built for modern wellness studios and professional environments, the vertical LedPro™ combines red and near-infrared light exposure, whole-body vibration training, and low-frequency sound programs in a compact upright format. Designed for comfortable standing use, it offers an immersive full-body wellness experience while saving space compared with traditional bed-style systems.

The vertical LedPro™ brings the signature red and near-infrared light experience of our bed-style LedPro™ into a freshly redesigned format with enhanced features for professional wellness settings.

Red and near-infrared light create a warm, targeted light environment designed to support relaxation, body care routines, and a refreshed feeling. The integrated vibration platform adds whole-body movement, helping users feel activated, energised, and more connected to their body during each session.

Low-frequency sound programs, inspired by Nogier frequencies, add a calming sensory layer to the experience. These resonant tones are designed to support relaxation and create a more immersive atmosphere.

Together, light, vibration, and sound create a multi-sensory wellness experience for users seeking relaxation, renewal, and a premium studio-based recovery ritual — all while standing comfortably.

CONTRAINDICATIONS

Whenever unsure about any listed or other potential contraindications, please make sure that the client consults their doctor or seeks medical advice from the treating medical professional before sessions.

Below are some of the possible contraindications:

- Pregnancy or suspected pregnancy.
- Pacemaker, implanted defibrillator, neurostimulator, insulin pump, or other active implanted electronic medical device, unless cleared by a healthcare professional and the device manufacturer.
- Recent surgery, recent fracture, open wounds, acute injury, severe sprain or strain, or recent trauma.
- Known or suspected blood clots, including deep vein thrombosis, pulmonary embolism, thrombophlebitis, or high clotting risk.
- Severe or unstable cardiovascular conditions, including recent heart attack, unstable heart disease, severe arrhythmia, uncontrolled high blood pressure, or any condition where exercise or vibration has been restricted.
- Photosensitivity disorder or any condition known to be triggered or worsened by light exposure.
- Use of photosensitising medication or topical products, unless cleared by a healthcare professional.
- Photosensitive epilepsy or history of light-triggered seizures, especially if using pulsed or flashing light modes.
- Active cancer, suspicious lesions, or use directly over a tumour area, unless specifically cleared by a qualified physician.
- Severe balance problems, fainting episodes, uncontrolled vertigo, severe dizziness, or inability to stand safely.
- Severe osteoporosis, unstable spinal conditions, recent disc injury, severe joint instability, or any condition where vibration or impact loading has been restricted.
- Recent eye surgery, serious retinal condition, or medical advice to avoid bright light or near-infrared light exposure.
- Active infection, fever, acute inflammation, unexplained swelling, or unexplained pain.
- Underage children.
- Any condition where a healthcare professional has advised avoiding vibration, exercise, bright light, heat exposure, low-frequency sound, or standing unsupported.

SAFETY PRECAUTIONS



WARNINGS AND PRECAUTIONS

INSTALLATION AND OPERATING AREA

Install the device only on a stable, hard, level, dry, and slip-resistant surface capable of supporting the device and user load.

Keep the area around the device clear so users can step on, step off, and move away safely.

Do not position the device where walls, furniture, accessories, cables, or other objects may restrict safe entry, exit, or operation.

Keep the power cord positioned so it cannot be tripped over, pinched, crushed, pulled, or damaged.

Do not use the device if it is unstable, rocking, visibly damaged, or not positioned securely.

Do not use the device outdoors.

Do not use the device near pools, showers, saunas, wet rooms, or other sources of water or excessive humidity.

ELECTRICAL SAFETY

Connect the device only to a mains supply that matches the voltage and frequency stated on the rating label.

Use only a properly grounded electrical outlet.

Do not remove, modify, or bypass the protective earth connection.

Ensure that the mains plug or disconnecting device remains accessible during use.

Do not use the device if the power cord, plug, socket, housing, control panel, lights, platform, or any electrical part appears damaged.

Do not operate the device with wet hands or while standing on a wet surface.

Do not expose the device, power cord, plug, or electrical components to water, liquids, excessive humidity, or cleaning fluids.

Do not immerse any part of the device in water or other liquid.

Do not use the device if it has been dropped, tipped over, damaged, exposed to water, or partially immersed in liquid.

If liquid enters the device, stop use immediately, disconnect the device from the mains supply, and have it inspected by qualified service personnel before further use.

Keep the power cord away from heated surfaces, sharp edges, moving parts, and walkways.

Do not pull, lift, carry, or move the device by the power cord.

Do not use extension cords, adapters, or multi-socket outlets unless they are suitably rated and permitted by the installation instructions.

Never insert objects into any opening, ventilation slot, gap, or moving area of the device.

If an object falls into the device, switch the device off, disconnect it from the mains supply, and contact qualified service personnel if the object cannot be safely removed.

If the device produces unusual noise, smell, smoke, excessive heat, sparks, or abnormal operation, stop use immediately, disconnect it from the mains supply, and contact qualified service personnel.

Do not open the housing or attempt to repair electrical components.

Electrical inspection and repair must be carried out only by qualified service personnel.

If the electrical outlet, circuit breaker, plug, or building supply appears faulty, contact a licensed electrician.

SAFETY PRECAUTIONS



WARNINGS AND PRECAUTIONS

MECHANICAL AND PLATFORM SAFETY

Do not exceed the stated maximum user weight or maximum load capacity.

Step onto and off the vibration platform carefully.

Use support handles or rails where provided.

Do not jump on the platform.

Do not use the platform for exercises or movements not described in the user instructions.

Do not sit, kneel, lie down, or place the head, chest, abdomen, or spine directly on the vibration platform unless the device is specifically designed and labelled for that use.

Never reach into, under, or around the vibration platform while the device is operating.

Do not tip, lift, move, or reposition the device while it is connected to power or in operation.

Do not place hands, feet, tools, accessories, or other objects under the device during operation.

Keep hands, feet, clothing, hair, jewellery, and loose items away from gaps, moving parts, ventilation openings, and underside areas.

Do not use the vibration platform if you cannot stand independently and safely.

Check regularly for loose fasteners, unusual vibration, abnormal noise, visible damage, worn parts, or instability.

Remove the device from service if damage, looseness, abnormal movement, or unsafe operation is observed.

USER POSITION AND FOOTWEAR

Maintain a balanced stance during vibration.

Keep knees slightly bent during vibration.

Do not lock the knees, elbows, or other joints during use.

Do not lean back heavily onto the heels during vibration.

Keep body weight controlled and avoid unstable body positions.

Use handles only for balance; do not hang from the handles.

Wear sports shoes. Do not use the vibration plate barefoot.

SAFETY PRECAUTIONS



WARNINGS AND PRECAUTIONS

GENERAL USER SAFETY

Consult a physician or qualified healthcare professional before use, especially if you have a medical condition, implanted device, injury, balance concern, light sensitivity, or are unsure whether the device is suitable for you.

This device is intended for general wellness use only and is not intended to diagnose, treat, cure, prevent, or alleviate any disease or medical condition.

Read the full user manual, including all warnings and safety instructions, before using the device.

The owner or operator must ensure that all users receive the necessary safety information before use.

Use the device only as described in the user instructions.

Do not use accessories, attachments, replacement parts, software, or modifications that are not approved by the manufacturer.

Stop use immediately if you experience pain, dizziness, nausea, shortness of breath, chest discomfort, headache, eye discomfort, skin irritation, excessive warmth, numbness, tingling, anxiety, or loss of balance.

Do not exceed the stated session time, intensity level, frequency of use, user weight, or load limits.

First-time users should begin with the lowest settings and shortest session duration.

Do not over-exert yourself or continue use to the point of exhaustion.

Allow sufficient rest between sessions.

Stay properly hydrated before and after use.

Do not use the device while tired, unwell, dizzy, dehydrated, or under the influence of alcohol, recreational drugs, or medication that may affect balance, alertness, or sensitivity.

Do not use immediately after a heavy meal, sauna, intense exercise, or other activity that may increase dizziness, dehydration, or discomfort.

Stop use immediately if you experience pain, dizziness, nausea, shortness of breath, chest discomfort, headache, eye discomfort, skin irritation, excessive warmth, numbness, tingling, anxiety, or loss of balance.

Retain the user manual for future reference.

Do not modify, alter, repair, bypass, disable, or tamper with any part of the device, including safety features, electrical components, controls, software, light modules, ventilation, guards, or the vibration platform, unless expressly authorised by the manufacturer.

USER SUPERVISION

This appliance is not to be used by persons, including children, with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning safe use of the appliance and understand the hazards involved.

Children must be supervised to ensure that they do not play with the appliance.

Children must not use, clean, or perform user maintenance on the appliance without appropriate supervision.

Keep children and pets away from the device during operation.

Commercial operators should supervise first-time users and ensure that users understand how to start, stop, and safely exit the device before each session.

Operators should screen users for contraindications and precautions before use.

Operators should not allow use by anyone who appears unable to stand safely, follow instructions, or respond appropriately during operation.

SAFETY PRECAUTIONS



WARNINGS AND PRECAUTIONS

VIBRATION PRECAUTIONS

Vibration may not be suitable for all users.

Users with balance, cardiovascular, neurological, musculoskeletal, spinal, implant-related, or injury-related concerns should consult a healthcare professional before use.

Start at the lowest vibration level and increase only if comfortable.

Stop the vibration function immediately if you feel dizziness, discomfort, pain, numbness, tingling, nausea, headache, or loss of balance.

Do not use the vibration function if a healthcare professional has advised you to avoid vibration, impact, or exercise.

Do not use the vibration function if you are unable to maintain a stable standing position.

Do not use dumbbells, weighted vests, resistance bands, or other added loads on the platform unless specifically permitted by your physician.

When additional load is used, reduce vibration intensity, session duration, and exercise complexity.

RED AND NEAR-INFRARED LIGHT PRECAUTIONS

Do not stare directly into the red or near-infrared light sources.

Use protective eyewear for red light and near-infrared applications. Ordinary sunglasses must not be relied upon.

Users with eye disease, recent eye surgery, retinal conditions, glaucoma, cataracts, macular degeneration, light sensitivity, or photosensitive conditions should consult a healthcare professional before use.

Do not use red/NIR light over skin that is injured, inflamed, infected, sunburned, recently treated, or unusually sensitive.

Use caution over tattoos, dark-pigmented areas, cosmetic products, or heat-sensitive skin.

Stop use if you feel uncomfortably warm, pain, irritation or unusual sensitiveness.

Users taking medication or using topical products that may increase light sensitivity should consult a healthcare professional before use.

Do not exceed the recommended session times.

Do not cover the light sources or block ventilation during operation.

SOUND PROGRAM PRECAUTIONS

Keep sound volume at a comfortable level.

Stop use if low-frequency sound causes dizziness, nausea, headache, anxiety, ear discomfort, or any unpleasant sensation.

Users with tinnitus, migraines, vestibular disorders, epilepsy, sound sensitivity, or neurological conditions should consult a healthcare professional before use.

OPERATING LIMITS

Do not operate the device continuously beyond the maximum session time specified by the manufacturer.

If the device becomes unusually hot, stop use, disconnect it from the mains supply, and allow it to cool before inspection.

Do not leave the device operating unattended.

Switch the device off or to standby mode after use.

SAFETY PRECAUTIONS



WARNINGS AND PRECAUTIONS

CLEANING AND USER MAINTENANCE

Switch off and disconnect the device from the mains supply before cleaning, inspection, troubleshooting, moving, or user maintenance.

After disconnecting, wait at least one minute before performing any permitted user maintenance.

Do not perform cleaning, inspection, or maintenance while the device is operating.

Clean the device only according to the cleaning instructions provided by the manufacturer.

Clean exterior surfaces only with a soft, slightly damp cloth unless otherwise instructed.

Do not spray cleaning solution directly onto the device, control panel, lights, ventilation openings, power cord, or plug.

Do not use sharp objects, abrasive pads, bristles, acid-based cleaners, solvents, or aggressive disinfectants unless approved by the manufacturer.

Clean electrical components, control panels, connectors, ventilation openings, and light modules only as instructed and never with excessive moisture.

Allow the device to dry completely before reconnecting it to the mains supply.

Do not perform service or maintenance beyond the user-maintenance tasks described in the instructions.

Servicing must be carried out only by qualified service personnel.

COMMERCIAL OPERATION

Operators must ensure that the device is installed, used, cleaned, inspected, and maintained according to the manufacturer's instructions.

Operators should ensure that only trained staff OR users operate, clean, move, or maintain the device where required.

Operators should keep the area around the device clean, dry, and free from trip hazards.


Operators should inspect the device regularly for damage, wear, loose parts, overheating, abnormal vibration, electrical faults, and unsafe operation.

Operators should remove the device from service if damage, abnormal operation, excessive vibration, overheating, or electrical faults are observed.


Operators should keep records of inspection, maintenance, service, and any reported incidents where appropriate.

UNBOXING INSTRUCTIONS

1.




2.




Unscrew and remove the panels.

3.




Ensure enough space for rolling the device off the pallet.

4.




Gently lift one end off the pallet.

5.




Gently lift the other end off the pallet.

6.



Lift up the device with caution. 3 people recommended.


7.



HEX 6

Unscrew the vibration plate transport screws on both sides. The vibration plate will drop to the floor.

8.



HEX 6

Retighten the screws on both sides.

INTRODUCTION TO FUNCTIONS

RED LIGHT AND NEAR-INFRARED

LedPro Vertical uses visible red light at 660 nm $\pm$ 10 nm and near-infrared light at 850 nm $\pm$ 10 nm to create a full-body upright light exposure experience. Red light creates the visible glow, while near-infrared light is mostly invisible and may be felt as gentle warmth depending on the selected intensity, session length and individual sensitivity.

The light intensity can be adjusted manually, allowing each session to feel softer or more intensive. Higher intensity is not automatically better: the overall experience depends on the combination of intensity, exposure time, skin sensitivity, and frequency of use. New users should start with lower settings and increase gradually.

In a wellness setting, red and near-infrared light exposure can support skin-focused routines, relaxation, post-activity rituals, and a refreshed feeling. It can be used alone for a calm light-focused session or combined with vibration and sound for a more immersive full-body experience.

VIBRATION TRAINING

LedPro Vertical includes an integrated vibration platform with adjustable frequency up to 40 Hz and amplitude up to 5 mm. Frequency controls how fast the platform moves, while amplitude describes the movement depth. In LedPro Vertical, amplitude is informational and adjusts automatically according to the selected vibration frequency, intensity level and weight of the user.

The vibration settings allow the session to range from gentle to more active. Lower settings are suitable for first-time users, warm-up style sessions, mobility routines, or softer wellness experiences. Higher settings create a stronger standing stimulus and should be increased gradually as the user becomes comfortable.

Vibration training can support body activation, balance engagement, muscle awareness, and active body-conditioning routines. Because the user is standing and responding to the platform movement, vibration sessions can make the experience more physically engaging and may increase energy expenditure compared with a fully passive session. When combined with a balanced lifestyle, this can support body-toning and fat-burning focused fitness routines.

When combined with red and near-infrared light, vibration turns a still light session into a more dynamic standing experience. The user remains comfortably in place while the body is engaged through posture, balance, and controlled platform movement.

NOGIER-INSPIRED SOUND FREQUENCIES

LedPro Vertical includes selectable sound frequencies inspired by the Nogier frequency concept: 73 Hz, 146 Hz, 292 Hz, 584 Hz, 1168 Hz, 2336 Hz, and 4672 Hz. These frequencies can be played through headphones as standalone tones or blended with music for a more atmospheric listening experience.

The sound function adds an audio layer to the session, helping shape the mood, rhythm, and energy of the experience. Users can choose calmer music for relaxation or more energetic music when they want to feel activated and ready to move. Standalone tones offer a focused listening experience, while music combinations can make the session feel smoother, more enjoyable, and more immersive.

Together, the three functions create a layered wellness ritual: light provides the visual and warming environment, vibration adds physical activation, and sound adds rhythm, focus, and atmosphere.


INSTRUCTIONS FOR USE

STARTUP PROCEDURE, STARTING AND STOPPING THE SESSION



1. Plug the device directly to the power outlet capable of supporting the electrical specifications. Do not use an extension wire.

2. Press the power on/ off button to turn on the device. Wait for the login screen to appear.




Power on/off

3. Tap on the user account. Then tap on the password field to reveal the keyboard and type in the password (the default password is "1") . Tap the **LOGIN** button to proceed.

4. For independent use, step inside the device and operate your session from the internal screen.

To initiate a session, choose a preferred **Session Setting** from one of the presets, or enter the settings manually.

5. Tap on the **START** icon to begin.



6. Tap on the **STOP** icon to end the session. When the session has ended, the device will revert to **STANDBY** mode.

7. The device can be left in **STANDBY** when not in use.

8. For a complete shutdown, tap **LOGOUT->SHUT-DOWN->YES**, then wait for the screen to turn off, and press the main power on/off button.

SESSION GUIDELINES

Before starting your LedPro Vertical session, choose clothing that feels comfortable and allows you to stand and move naturally. For the most direct red and near-infrared light exposure, wear clothing that leaves the target areas of the skin uncovered. Depending on your comfort and privacy, suitable options may include swimwear, sports underwear, shorts, a sports bra, or light activewear.

Wear clean, stable trainers, running shoes, or similar supportive footwear when using the vibration platform. Avoid high heels, unstable shoes, wet footwear, slippery socks, or anything that may affect balance. Remove jewellery, watches, loose accessories, and items from pockets before the session. Tie back long hair and make sure clothing, straps, and accessories cannot contact the vibration platform, ventilation openings, or any moving areas.

Avoid applying heavy oils, lotions, reflective creams, or products that may make the skin or platform slippery before use. Clean, dry skin is recommended for comfort, hygiene, and consistent light exposure.

Use the supplied or approved eye protection during the session and do not look directly into the red or near-infrared light sources for an extended time. In addition to eye protection, we recommend keeping your eyes closed when possible. Headphones may be used for the sound programs, either with stand-alone frequencies or combined with music. Choose calmer audio for a relaxing session or more energetic music for a more active, body-engaging experience.

The recommended maximum session duration is 8-12 minutes, 2-3 times a week. More frequent or longer use is not necessarily better, as the overall experience depends on light intensity, vibration level, session duration, and individual sensitivity.

SESSION PRESETS



LedPro Vertical includes pre-programmed session settings designed to offer a simple starting point for different wellness and fitness goals. Each preset combines the three main functions of the device: red and near-infrared light exposure, whole-body vibration training, and Nogier-inspired sound frequencies.

The presets are built around the principle that each function can shape the session in a different way. Light intensity controls the strength of the red and near-infrared light exposure. Vibration frequency and intensity define how gentle or active the standing experience feels. Sound frequency adds an audio layer that can support the overall mood, rhythm, and focus of the session.


Each setting uses a different balance of these elements. Some presets place more emphasis on light exposure, while others create a more active vibration-based experience or a stronger sensory atmosphere through sound. The presets are intended as guided starting points and can be adjusted manually depending on user comfort, experience level, and session preference.

New users should begin with lower manual settings or a comfortable preset, then increase intensity gradually only if the session feels safe and comfortable. Sessions should always be used according to the contraindications, precautions, and safety instructions provided in this manual.


LedPro Vertical Predetermined Session Settings

Session Setting	Time (min)	Vibration Frequency (Hz)	Vibration Intensity (%)	Light Intensity (%)	Sound Frequency (Hz)
Body Toning	10	9	100	75	1168
Skin Toning	10	4	80	80	292
Inflammation Reduction	10	30	75	70	292
Mobility	10	16	80	100	4672
Stress Relief	10	29	100	40	2336

MANUAL SESSION SETTINGS




Tap on the individual function icons to set session parameters.



1. Session time

Adjust the slider to change the duration of the session. We recommend following the max recommended session duration set out in [Session Guidelines](#).





2. Vibration settings

Choose between a simplified and an advanced vibration settings menu.

The simplified menu highlights the focus area which the vibration plate will affect most along with an intensity adjustment.

The advanced menu offers a numeric approach to vibration training frequency and intensity.

Plate Synchronisation determines whether both sides of the platform move in the same phase or in opposing phases. In Synchronized mode, both sides rise and lower together. In Alternating mode, one side rises while the other lowers, producing a side-to-side vibration pattern.



3. Light intensity

Adjust the slider to change the light intensity.

4. Audio setting

Choose between the preferred sound frequency and it's volume. Additionally, music can be played at the same time.



MAINTENANCE

The LedPro™ device does not require any particular maintenance besides a basic cleaning routine. Turn off the device before any cleaning procedures and unplug it. Disinfect the stability handles after use. Wipe off the device itself with a damp soft microfibre cloth if needed. For wiping the doors, we recommend alcohol-free disinfection and cleaning wipes that are suitable for sensitive surfaces (for example, Sterisept Wipes by Chemi-Pharm. Alternatively, use a soft damp microfibre towel. If the device appears to function abnormally, do not attempt to open the device and repair it. Please contact CTN support at service@ctn.fi.

WARRANTY INFORMATION

The standard limited warranty for this product is 24 months. For additional information, please refer to the CTN limited warranty document.

WORLD'S #1 IN SPORTS RECOVERY



CTN GROUP - IRELAND - FINLAND - ESTONIA - SPAIN

**CRYOTECH FINLAND OY - JAAKONKATU 2 - 01620 VANTAA - FINLAND
CRYOTECH NORDIC OÜ - HÄRGMÄE 20 - 13525 TALLINN - ESTONIA
CTN RECOVERY STUDIO - PASAJE MUIXARA 11 - 03530 ALICANTE - SPAIN**

SALES@CTN.FI ADMIN@CTN.FI SERVICE@CTN.FI