

# X°CRYO™ APPLICATOR GUIDE



## Cryo°Sleeve Ø150 M Ø200 L

Designed for the hands, knees and elbows. Available in M & L sizes.



## Cryo°Massager Ø35 S Ø70 M

The favourite of many physiotherapists and masseurs, for providing a simultaneous cooling and massage.



## Cryo°Spot Ø12 S Ø20 M

Suitable for small joint areas or skin lesions. Standard two sizes.



## Cryo°Foot

Designed for ankles, soles and Achilles tendons.



## Cryo°Mask

Refresh your complexion with cryofacials, which tighten pores, reduce puffiness and promote collagen production.


\*DIAMETERS Ø IN MILLIMETERS

# CTN.FI



## Cryo°Meter®

Designed to connect to the X°CRYO™ hose for temperature monitoring and ease of use. Compatible with the Cryo°Cups and the Cryo°Nozzle.



## Cryo°Nozzle Ø25 L

The Cryo°Nozzle helps target smaller areas and is compatible with the Cryo°Meter™.



## Cryo°Cup Ø60 XS Ø90 S Ø130 M Ø180 L

The Cryo°Cup applicators are suitable for a variety of areas, compatible with the Cryo°Meter™.

